

Use and Disposal of single-patient Blood Glucose Monitoring Kit for C-diff patient/resident

- *The assigned blood glucose meter will be left in the patient's room and used only for him/her during the time they have C-diff.*
- *Do NOT remove this device from the C-diff resident's room*
- *Once the patient/resident has been cleared from C-diff isolation, you must dispose of the blood glucose meter.*



1



Wash hands thoroughly with soap and water. Wear a clean pair of gloves.

2



Make sure meter is off and a test strip is not inserted.

5



Allow meter to air dry.

3



With approved wipe, rub the entire outside of the meter using 3 circular wiping motions with moderate pressure on the front, back, left side, right side, top and bottom of the meter. Repeat as needed until all surfaces are visibly clean. Discard used wipe.

6



Dispose of used wipes and gloves. Wash hands thoroughly with soap and water.

4



Using a fresh wipe, make sure that all outside surfaces of the meter remain wet for 2 minutes. **Make sure no liquids enter the test port or any other opening.**



The following product has been approved to clean the meter:

- Super Sani-Cloth Wipes (EPA Reg. No. 9480-4), manufacturer: PDI Healthcare.
- The meter has been validated for 5,200 cleaning and disinfection cycles.

Tyson Bioresearch, Inc.

5F., No. 16, 18, 20, 22, Kedong 3rd Rd.,
Zhunan Township, Miaoli County 35053,
Taiwan
+886-37-585998 www.tysonbio.com

Customer Service

U.S.A. +1-800-948-3894
24 hours a day, 7 days a week
Taiwan +886-37-585998
GMT+8, 8:30 ~ 17:30, Monday ~ Friday