

Cleaning and Disinfecting Guide

A GUIDE FOR MULTIPLE PATIENT USE FACILITIES

The meter should be disinfected after use on each patient. This Blood Glucose Monitoring System may only be used for testing multiple patients when Standard Precautions and the manufacturer's disinfection procedures are followed.

All parts of the glucose monitoring system should be considered potentially infectious and are capable of transmitting bloodborne pathogens between patients and healthcare professionals. For more information, refer to "Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings 2007", <http://www.cdc.gov/hicpac/2007ip/2007isolationprecautions.html>



1



Wash hands thoroughly with soap and water. Wear a clean pair of gloves.

2



Make sure meter is off and a test strip is not inserted.

5



Allow meter to air dry.

3



With approved wipe, rub the entire outside of the meter using 3 circular wiping motions with moderate pressure on the front, back, left side, right side, top and bottom of the meter. Repeat as needed until all surfaces are visibly clean. Discard used wipe.

6



Dispose of used wipes and gloves. Wash hands thoroughly with soap and water.

4



Using a fresh wipe, make sure that all outside surfaces of the meter remain wet for 2 minutes. **Make sure no liquids enter the test port or any other opening.**



The following product has been approved to clean the meter:

- Super Sani-Cloth Wipes (EPA Reg. No. 9480-4), manufacturer: PDI Healthcare.
- The meter has been validated for 5,200 cleaning and disinfection cycles.

Tyson Bioresearch, Inc.

5F., No. 16, 18, 20, 22, Kedong 3rd Rd.,
Zhunan Township, Miaoli County 35053,
Taiwan
+886-37-585998 www.tysonbio.com

Customer Service

U.S.A. +1-800-948-3894
24 hours a day, 7 days a week
Taiwan +886-37-585998
GMT+8, 8:30 ~ 17:30, Monday ~ Friday