

Waterproof Foam Dressing with Adhesive Border

Waterproof Foam Dressing with Adhesive Border is a moderately absorbent, waterproof sterile foam dressing. It is structured to retain moderate exudate, while maintaining a moist wound environment.

Waterproof Foam Dressing with Adhesive Border has a waterproof adhesive backing and a sterile foam cushioning that easily adheres to the skin and seals the wound bed, protecting against bacterial contamination. The waterproof outer layer of the dressing prevents water or other liquids from penetrating the wound site, making it ideal for use in wet environments. The adhesive border provides a secure and comfortable fit, helping to keep the dressing in place and preventing it from slipping or rolling.



ENHANCED PATIENT COMFORT

Foam dressings are designed to be highly absorbent, while still being gentle to the skin, with minimal skin stripping upon removal.



Not made with natural rubber latex.



Indications:

A sterile primary dressing that may be used for the management of moderately to highly exuding wounds, such as leg and foot ulcers, pressure ulcers, diabetic foot ulcers, skin abrasions, surgical wounds, donor sites, 1st and 2nd degree burns.

Features & Benefits:

- Provides a balanced moisture environment to promote autolytic debridement
- Minimizes maceration and promotes healing with excellent exudate absorption
- Maintains a moist environment for optimal wound healing
- Conforms easily to awkward-to-dress areas
- Not manufactured with natural rubber latex
- Waterproof for showers or baths

Available Sizes:

Products	Size	Reorder #	Units
Foam Dressing/Border	4in x 4in	PCHF3005	10/Box, 40 Boxes/Case
Foam Dressing/Border	6in x 6in	PCHF3006	10/Box, 20 Boxes/Case

Compare to ALLEVYN ADHESIVE • TIELLE ADHESIVE • BIATAIN • MEPILEX TRANSFER / BORDER • OPTIFOAM ADHESIVE • PERMAFOAM, DERMA LEVIN