

Soft Silicone Foam Dressing without Border

Soft Silicone Foam Dressing without Border is for the management of moderate to heavily exuding wounds and consists of a flexible, highly absorbent, advanced soft silicone foam layer. The absence of a border allows the dressing to be easily cut and customized to fit the specific wound size and shape. The foam pad consists of a vapor-permeable, moisture-proof outer film, and a silicone adhesive that adheres to the surrounding skin, allowing exudate to be absorbed. The dressing is easy to apply and remove minimizing pain, maceration and the risk of secondary damage to the wound area.



ENHANCED PATIENT COMFORT

The superb absorption of foam and the gentle adhesion of silicone minimizes pain and trauma during dressing change.

Not made with natural rubber latex.



Available Sizes:

Products	Size	Reorder #	Units
Soft Foam Dressing	2in x 2in	PCHF3100	10/Box, 40 Boxes/Case
Soft Foam Dressing	3in x 3in	PCHF3101	10/Box, 40 Boxes/Case
Soft Foam Dressing	4in x 4in	PCHF3102	10/Box, 40 Boxes/Case
Soft Foam Dressing	4in x 8in	PCHF3103	10/Box, 20 Boxes/Case
Soft Foam Dressing	6in x 6in	PCHF3104	10/Box, 20 Boxes/Case
Soft Foam Dressing	6in x 8in	PCHF3105	5/Box, 40 Boxes/Case
Soft Foam Dressing	8in x 8in	PCHF3106	5/Box, 40 Boxes/Case

Indications:

A primary dressing that may be used for the management of moderate to highly exuding wounds, such as leg and foot ulcers, pressure ulcers, traumatic wounds (e.g., skin tears and secondary healing wounds), under compression, graft and donor sites.

Features & Benefits:

- Safe, easy to apply and remove, and can be sized to fit the wound
- Highly absorbent for moderate to heavy exudate
- Maintains a moist environment for wounds and promotes wound healing
- Seals and prevents the introduction of foreign particles and debris into the wound area
- Effective and sustained protection for several days, allowing for optimal healing without disturbance

Compare to MEPILEX • MEPILEX BORDER
MEPILEX BORDER LITE • ALLEVYN GENTLE BORDER
ALLEVYN GENTLE BORDER LITE • OPTIFOAM
GENTLE BORDER • AQUACEL FOAM • BIATAIN
TIELLE ESSENTIAL • KERRAFOAM • POLYDERM GTL
COMFORTFOAM • COMFORTFOAM BORDER
COMFORTFOAM BORDER LITE

