PROCURE HEALING

Soft Silicone Foam Dressing-LITE

Slimmer Version with Adhesive Border

Soft Silicone Foam Dressing-LITE with Border is a slimmer foam dressing for the management of low to moderate exuding wounds and consists of a thin, absorbent soft silicone foam layer. This highly comfortable and flexible foam pad consists of a vapor-permeable, moisture-proof outer film, and a silicone adhesive with a low profile border that adheres to the surrounding skin, allowing exudate to be absorbed, while preventing leakage. The dressing is easy to apply and remove minimizing pain, maceration and the risk of secondary damage to the wound area.



Products	Size	Reorder #	Units
Foam Dressing LITE/Border	2in x 2in	PCHF3200	10/Box, 80 Boxes/Case
Foam Dressing LITE/Border	3in x 3in	PCHF3201	10/Box, 80 Boxes/Case
Foam Dressing LITE/Border	4in x 4in	PCHF3202	10/Box, 40 Boxes/Case



HEALING

10_{PCS}

PROCURE Soft Silicone Foam Dressing-LITE

PR

Soft Silicone Foam Dressing-LITE

Waterproof and bacterial resistant
Low profile border for added adhesion
and flexibility
Helps maintain a moist environment
to promote healing

4 in X 4 in (10 x 10cm) 10 Dressings / PAD: 2.5 in x 2.5 in (6.3 x 6.3 cm)

Indications:

A primary dressing that may be used for the management of mild to moderate exuding wounds, especially for the hard to dress area, such as leg and foot ulcers, pressure ulcers, traumatic wounds, partial thickness burns.

Features & Benefits:

- Thin, comfortable, safe, easy to apply and remove
- Sterile, self-adhesive and showerproof dressing with a low profile border to seal the wound area and prevent exudate leakage
- Superior flexibility and comfortable adherence to skin, while minimizing pain, trauma and wound maceration
- Effective and gentle wound protection for several days

Compare to MEPILEX • MEPILEX BORDER MEPILEX BORDER LITE • ALLEVYN GENTLE BORDER ALLEVYN GENTLE BORDER LITE • OPTIFOAM GENTLE BORDER • AQUACEL FOAM • BIATAIN TIELLE ESSENTIAL • KERRAFOAM • POLYDERM GTL COMFORTFOAM • COMFORTFOAM BORDER COMFORTFOAM BORDER LITE