

Advanced Wound Care

Soft Silicone Foam Dressing

with Adhesive Border

Soft Silicone Foam Dressing with Adhesive Border is for the management of moderate to heavily exuding wounds and consists of a flexible, highly absorbent, advanced soft silicone foam layer. The foam pad consists of a vapor-permeable, moisture-proof outer film, and a silicone adhesive border that adheres to the surrounding skin, allowing exudate to be absorbed, while preventing leakage. The dressing is easy to apply and remove minimizing pain, maceration and the risk of secondary damage to the wound area.

Features & Benefits:

- Sterile, self-adhesive and showerproof dressing with a low profile border to seal the wound area
- Highly absorbent for moderate to heavy exudate, while still allowing breathability to promote healing
- Supports and accelerates moist wound healing
- Conforms and protects fragile tissue and is suitable for wounds on various body parts
- · Provides thermal insulation, absorbing high levels of drainage
- Effective and sustained protection for several days, allowing for optimal healing without disturbance

Available Sizes:

| Products | Size | Reorder# | Units |
|------------------|--------------------------|----------|-----------------------|
| Soft Foam/Border | 2in x 2in | PCHF3300 | 10/Box, 80 Boxes/Case |
| Soft Foam/Border | 2in x 5in | PCHF3301 | 10/Box, 40 Boxes/Case |
| Soft Foam/Border | 3in x 3in | PCHF3302 | 10/Box, 80 Boxes/Case |
| Soft Foam/Border | 4in x 4in | PCHF3303 | 10/Box, 40 Boxes/Case |
| Soft Foam/Border | 4in x 8in | PCHF3304 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | 4in x 12in | PCHF3305 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | 6in x 6in | PCHF3306 | 10/Box, 20 Boxes/Case |
| Soft Foam/Border | 6in x 8in | PCHF3307 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | 7.2in x 7.2in | PCHF3308 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | Heel-5.3in x 6in | PCHF3309 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | Sacral(sm)-7.2in x 7.2in | PCHF3310 | 5/Box, 40 Boxes/Case |
| Soft Foam/Border | Sacral(lg)-9in x 9in | PCHF3311 | 5/Box, 40 Boxes/Case |



Indications:

- Exuding and Traumatic Wounds (Effectively Seals Cavity Wounds)
- Pressure and Diabetic Ulcers
- Venous and Arterial Ulcers
- Partial and Full Thickness Wounds (including Second and Third Degree Burns)

ENHANCED PATIENT COMFORT

The superb absorption of foam and the gentle adhesion of silicone minimizes pain and trauma during dressing change.





Not made with natural rubber latex.

Compare to MEPILEX • MEPILEX BORDER • MEPILEX BORDER LITE • ALLEVYN GENTLE BORDER • ALLEVYN GENTLE BORDER LITE • OPTIFOAM GENTLE BORDER • AQUACEL FOAM • BIATAIN • TIELLE ESSENTIAL • KERAFOAM • POLYDERM GTL • COMFORTFOAM BORDER • COMFORTFOAM BORDER • COMFORTFOAM BORDER LITE

DEPTH/STAGE PARTIAL/FULL

1 2 3 4

