

# **Silver Silicone Foam Dressing**

without Border

Silver Silicone Foam Dressing without Border features a silver silicone island dressing that provides a vapor permeable waterproof backing for long-term protection. The absence of a border allows the dressing to be easily cut and customized to fit the specific wound size and shape. The gentle adhesion of silicone especially minimizes pain and trauma when dressing changing.

Silver Silicone Foam Dressing without Border contains 0.25–0.35 mg/cm2 silver ions that promote effective antibacterial protection. This dressing absorbs wound exudate and releases silver ions within the dressing in the presence of wound fluid to help reduce bacterial colonization of the dressing. It also assists in the maintaining a moist environment for optimal wound healing.

# ENHANCED PATIENT COMFORT

The superb absorption of foam and the gentle adhesion of silicone minimizes pain and trauma during dressing change.

Not made with natural rubber latex.

Ag

#### **Available Sizes:**

Products	Size	Reorder #	Units
AG Foam Dressing	2in x 2in	PCHF3400	10/Box, 40 Boxes/Case
AG Foam Dressing	4in x 5in	PCHF3401	10/Box, 40 Boxes/Case
AG Foam Dressing	6in x 6in	PCHF3402	10/Box, 20 Boxes/Case



## **Indications:**

A primary dressing that may be used for the management of moderately to highly exuding wounds, such as leg and foot ulcers, pressure ulcers, diabetic foot ulcers, skin abrasions, surgical wounds, donor sites, 1st and 2nd degree burns.

## **Features & Benefits:**

- · Absorbs and retains exudate
- Allows for observation and wound drainage assessment without disturbance during the healing process
- · Sterile and self-adhesive
- Effective long-term use
- Protects fragile tissue and maintains a moist and soft wound environment





Compare to Mepilex AG, Allevyn AG, Biatain AG, Optifoam Gentle AG, Comfort Foam AG