

Silver Silicone Foam Dressing

with Adhesive Border

Silver Silicone Foam Dressing with Adhesive Border features a silver silicone island dressing that provides a vapor-permeable waterproof backing for long-term protection. The gentle adhesion of silicone especially minimizes pain and trauma when dressing changing.

Silver Silicone Foam Dressing with Adhesive Border contains 0.25 – 0.35 mg/cm2 silver ions that promote effective antibacterial protection. This dressing absorbs wound exudate and releases silver ions within the dressing in the presence of wound fluid to help reduce bacterial colonization of the dressing. It also assists in the maintaining a moist environment for optimal wound healing.

ENHANCED PATIENT COMFORT

The superb absorption of foam and the gentle adhesion of silicone minimizes pain and trauma during dressing change.

Not made with natural rubber latex.

Ag

Available Sizes:

Products	Size	Reorder#	Units
AG Foam Dressing/Border	2in x 2in	PCHF3500	10/Box, 80 Boxes/Case
AG Foam Dressing/Border	4in x 4in	PCHF3501	10/Box, 40 Boxes/Case
AG Foam Dressing/Border	6in x 6in	PCHF3502	10/Box, 20 Boxes/Case

Advanced Wound Care PROCURE HEALING Silver Silicone Foam Dressing with Adhesive Border Antimicrobial silver effect Highly absorbent and waterproof Moisture wound healing Ag Sin x Sin (15x15cm) Moreometric Highly absorbent and waterproof Month Adhesive Border Ag Sin x Sin (15x15cm) Moreometric Highly Adhesive Ag Recompting High 43in x 3in (11x11cm) Ag Sin x Sin (15x15cm) Moreometric Highly Adhesive Ag Recompting High 43in x 3in (11x11cm) Moreometric Highly Adhesive Ag Ag Sin x Sin (15x15cm) Moreometric Highly Adhesive Ag Moreometric Highly Ag Moreometri

Indications:

A primary dressing that may be used for the management of moderately to highly exuding wounds, such as leg and foot ulcers, pressure ulcers, diabetic foot ulcers, skin abrasions, surgical wounds, donor sites, 1st and 2nd degree burns.

Features & Benefits:

- Seals the wound area while absorbing and retaining exudate
- Allows for observation and wound drainage assessment without disturbance during the healing process
- Sterile, self-adhesive, and showerproof
- Flexibility and adherence to skin
- Effective long-term use for up to 7 days
- Protects fragile tissue and maintains a moist and soft wound environment

DEPTH/STAGE PARTIAL/FULL

1 2 3 4



Compare to Mepilex AG, Allevyn AG, Biatain AG, Optifoam Gentle AG, Comfort Foam AG