

Non-Skid **Shoe Covers**

Procure Shoe Covers are protective coverings designed to be worn over shoes, providing a barrier between footwear and the environment. Made of Non-Woven spunbonded polypropylene, they are light and breathable yet fluid-resistant. Commonly used in a variety of settings including healthcare, construction, cleanrooms and other environments where maintaining cleanliness and preventing contamination is crucial.

Shoe covers help minimize the transfer of dirt, dust, and other substances from shoes. They are disposable, convenient and easy to slip on and off, and a practical solution for maintaining safety standards. The Non-skid material helps with ease of use with movement and traction.



Features:

- Made of fluid-resistant and breathable Non-Woven spunbonded polypropylene material to prevent fluid penetration and increased comfort.
- Disposable and single-use to reduce the risk of cross-contamination and infection transmission.
- Non-skid for traction
- Not made with natural rubber latex (Latex-Free).



Product	Reorder #	Units
Shoe Covers	PC631	100 Per Pack(50 pairs)/Box, 10 Packs/Case

Visit Us Online at:
www.ProCureProducts.com