



Tuberculin "TB" Safety Syringes *Sterile*

Tuberculin Safety Syringes are used for administering small volumes of medication, typically for tuberculin skin tests or other diagnostic purposes. The Safety syringe contains a protective sheath that covers the needle before and after use to help prevent needlestick injuries, on completion of the injection, the safety mechanism locks the needle in the barrel. These syringes have a small capacity, usually ranging from 0.5 mL to 1 mL, and are marked with fine graduations to enable precise measurement of medication. The fine gauge needle provides a smoother insertion for an improved comfortable injection experience.



Features & Benefits:

- Not made with natural rubber latex
- Sterile - individually packaged
- Easy to secure pull-and-twist safety mechanism
- Pre-attached ultra-sharp stainless steel needle
- Transparent barrel with large-font graduations allows for accurate measurement
- Low dead-space design

Visit Us Online at:
www.ProCureProducts.com

Product	Size	Needle	Reorder #	Units
Tuberculin Syringe	1ml	25g x 5/8"	PCSS730	100 Syringes/Box, 8 Boxes/Case
Tuberculin Syringe	1ml	27g x 1/2"	PCSS735	100 Syringes/Box, 8 Boxes/Case
Tuberculin Syringe	1ml	28g x 1/2"	PCSS740	100 Syringes/Box, 8 Boxes/Case