

### Instant Cold Packs

**Procure Instant Cold Packs** are for use whenever instant cold therapy is required. Immediate, colder-than-ice treatment for relief of pain and swelling caused by sprains, strains, contusions, toothaches, insect bites, minor injuries and minor burns. No refrigeration – squeeze and shake to activate. Once activated, these packs will consistently maintain a cold temperature for a full 20 minutes of clinically recommended application. Dispose of the ice pack after use, single use



#### Benefits:

- Instant relief minor injuries
- Relief minor pain and swelling
- Helps with minor burns
- Disposable, single use only
- 6x9 inches




Visit Us Online at:  
[www.ProCureProducts.com](http://www.ProCureProducts.com)

contusions, minor burns, toothaches and insect bites.

Instructions for Use:



Grasp outer edges of pack at the arrows above. SQUEEZE firmly to rupture inner pouch.



SHAKE pack well to mix contents after activation.



APPLY to patient.



DISPOSE of used pack in proper receptacle.

Product	Size	Reorder #	Units
Cold Packs	6" x 9"	PC1210	24 Cold Packs/Case