

## Hydrocolloid-THIN *Sterile Wound Dressing with Border*

**ProCure Hydrocolloid Wound Dressings** are a type of wound care product designed to promote healing in a moist environment while providing a protective seal around the wound. Hydrocolloid dressings are composed of gel-forming agents, such as gelatin or carboxymethylcellulose, which absorb helps to maintain a moist environment conducive to healing and facilitates autolytic debridement, where the body's own enzymes break down dead tissue.

The border is made of an adhesive material that sticks to the surrounding skin which helps to protect the wound from external contaminants and friction, reducing the risk of infection and promoting comfort during wear.



### Indications:

Hydrocolloid dressings with borders are particularly suitable for low to moderately exuding wounds, such as pressure ulcers, minor burns, and superficial wounds.

### ENHANCED PATIENT COMFORT

*Designed to allow the wound to breathe, minimizing skin and tissue maceration while maintaining a moist environment to encourage healing.*

*Not made with natural rubber latex.*



### Features & Benefits:

- Provides thermal insulation and protection
- Comfortable and conformable
- Helps maintain a moist environment to promote healing
- Absorbent
- Long-lasting and pain reducing
- Available in various shapes and sizes
- Sterile, single use device

### Available Sizes:

Products	Size	Reorder #	Units
Hydrocolloid-THIN Dressing	2in x 2in	PCHH3102	20/Box, 20 Boxes/Case
Hydrocolloid-THIN Dressing	4in x 4in	PCHH3104	10/Box, 40 Boxes/Case
Hydrocolloid-THIN Dressing	6in x 6in	PCHH3106	5/Box, 20 Boxes/Case
Hydrocolloid-THIN, SACRAL	6in x 7in	PCHH3107	5/Box, 20 Boxes/Case
Hydrocolloid-THIN, HEEL/ELBOW	3.5in x 5in	PCHH3355	10/Box, 20 Boxes/Case

Visit Us Online at:  
[www.ProCureProducts.com](http://www.ProCureProducts.com)

Compare to DERMAFILM • COMFEEL • ULTEC •  
DUODERM • EXUDERM • NU-DERM