

PREMIUM Non-Woven Sponges

Premium Non-Woven Gauze Sponges are a high quality rayon/ polyester blend that offers the perfect balance of strength and absorbency. Enhanced softness, breathability and comfort.

The nonwoven design keeps the sponges lint-free, lowering the risk of fibers remaining in the wound, which can be critical for preventing infection and promoting efficient healing. They provide superior absorbency and fast wicking action ideally suited for a variety of medical uses.



Features & Benefits:

- Rayon and polyester blend.
- Available in sterile and non-sterile options
- Virtually non-adherent and lint-free.
- Faster wicking ability promotes a healthy healing environment.

Visit Us Online at:
www.ProCureProducts.com

Product	Size	Reorder #	Units
Non-Woven Sponges - 4PLY, Sterile	2in x 2in	PCNG224S	2 Sponges/Pouch, 25 Pouches/Box, 30 Boxes/Case
Non-Woven Sponges - 4PLY, Non-Sterile	2in x 2in	PCNG224N	200 Sponges/Bag, 20 Bags/Case
Non-Woven Sponges - 4PLY, Sterile	3in x 3in	PCNG334S	2 Sponges/Pouch, 40 Pouches/Box, 30 Boxes/Case
Non-Woven Sponges - 4PLY, Non-Sterile	3in x 3in	PCNG334N	200 Sponges/Bags, 10 Bags/Case
Non-Woven Sponges - 4PLY, Sterile	4in x 4in	PCNG444S	2 Sponges/Pouch, 50 Pouches/Box, 12 Boxes/Case
Non-Woven Sponges - 4PLY, Non-Sterile	4in x 4in	PCNG444N	200 Sponges/Bags, 10 Bags/Case