

Gelling Fiber Dressing *Sterile CMC Wound Dressing*

ProCure Gelling Fiber Dressing, Sterile CMC Wound Dressing is a highly absorbent wound care solution made from CMC (carboxymethylcellulose). It is designed to gel on contact with exudate, creating a moist wound environment that promotes healing while managing fluid levels and protecting the wound from infection.



ENHANCED PATIENT COMFORT

An ultra-absorbent dressing that transforms into a strong clear gel upon contact with wound exudate. This gel captures debris and bacteria, securely holding them within the fiber dressing.

Not made with natural rubber latex.



Features & Benefits:

- Sterile CMC wound dressing for enhanced fluid control
- Gels upon contact with exudate, forming a moist healing environment
- Ideal for moderate to heavily exuding wounds
- Promotes faster healing by maintaining proper moisture levels
- Non-adherent and easy to remove without disrupting the wound

Available Sizes:

Products	Size	Reorder #	Units
Gelling Fiber CMC Dressing	2in x 2in	PCHCMC22	1 Roll/Pouch, 10 Pouches/Box, 40 Boxes/Case
Gelling Fiber CMC Dressing	4in x 5in	PCHCMC45	1 Roll/Pouch, 10 Pouches/Box, 40 Boxes/Case
Gelling Fiber CMC Dressing	6in x 6in	PCHCMC66	1 Roll/Pouch, 5 Pouches/Box, 40 Boxes/Case
Gelling Fiber CMC Dressing -ROPE	.75in x 18in	PCHCMC075	1 Roll/Pouch, 5 Pouches/Box, 40 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com