

Combine ABD Pads

Procure Combine ABD pads are large and very useful for wounds on the abdomen. A combine ABD pad is for heavily bleeding or exuding wounds. The soft, non-woven outer layer works to wick moisture and fluids away from the wound and into the inner layer. The sealed ends help prevent the bandage from seeping.

They can be either a primary or secondary dressing. Since combine pads do not have adhesive, they need to be secured by a tape or wrap.



Features & Benefits:

- Super absorbent abdominal pads
- Great for wound dressing and wound packing
- Outer layer is a soft, non-woven material
- Fluffy inner material absorbs and disperses fluids
- Sealed edges prevent leakage
- Available in sterile and non-sterile options
- Not made with natural rubber latex (Latex-Free)

Visit Us Online at:
www.ProCureProducts.com



Product	Size	Reorder #	Units
Combine ABD Pads - Sterile	5" x 9"	PCABD59S	1 Pad/Pouch, 25 Pouches/Box, 16 Boxes/Case
Combine ABD Pads	5" x 9"	PCABD59N	25 Pads/Bag, 16 Bags/Case
Combine ABD Pads - Sterile	8" x 10"	PCABD810S	1 Pad/Pouch, 20 Pouches/Box, 16 Boxes/Case
Combine ABD Pads	8" x 10"	PCABD810N	25 Pads/Bag, 16 Bags/Case