

Transparent Film Dressing

Sterile Semipermeable Adhesive

ProCure Transparent Film Dressing is a sterile, semipermeable adhesive dressing with a framed border for easy application and secure adhesion. Designed to protect against water, bacteria, and contaminants while allowing moisture vapor exchange, it supports optimal healing conditions. The transparent film enables continuous wound observation without removal, and the integrated label provides space for documenting date and other important information—ideal for clinical efficiency and wound care management.



ENHANCED PATIENT COMFORT

This non-absorbent, moisture-vapor permeable transparent film dressing features a picture frame border for easy application. It is ideal for non-infected wounds with minimal drainage or as a secondary dressing. The dressing helps maintain a moist wound environment and prevents bacterial contamination.



Not made with
natural rubber latex.

Features & Benefits:

- Transparent film for easy wound monitoring
- Semipermeable to allow moisture vapor and oxygen exchange
- Sterile and adhesive for secure placement
- Ideal for cuts, abrasions, and post-operative wounds
- Waterproof barrier to protect against contaminants

Visit Us Online at:
www.ProCureProducts.com

Available Sizes:

Products	Size	Reorder #	Units
Transparent Film Dressing	2.375in x 2.75in	PCHTD23	1/Pouch, 100 Pouches/Box, 40 Boxes/Case
Transparent Film Dressing	4in x 5in	PCHTD45	1/Pouch, 50 Pouches/Box, 40 Boxes/Case
Transparent Film Dressing	4in x 10in	PCHTD410	1/Pouch, 10 Pouches/Box, 40 Boxes/Case
Transparent Film Dressing	6in x 8in	PCHTD68	1/Pouch, 10 Pouches/Box, 40 Boxes/Case
Transparent Film Dressing	8in x 12in	PCHTD812	1/Pouch, 10 Pouches/Box, 40 Boxes/Case