

Transparent Film with Island Dressing

Sterile Semipermeable Adhesive Dressing with Absorbent Pad

ProCure Transparent Film with Island Dressing is a sterile semipermeable adhesive with a framed border that combines a breathable film and an absorbent island pad. It protects wounds while allowing moisture vapor exchange. The film serves as a barrier to water, bacteria, and contaminants, while the island pad absorbs exudate and cushions the wound. The framed adhesive border ensures secure, easy application and reduces edge lift, making it ideal for post-operative sites, abrasions, and low to moderately exuding wounds.



**ENHANCED
PATIENT COMFORT**
This transparent dressing is designed to allow moisture vapor to escape while helping to prevent bacterial contamination. It easily conforms to difficult areas and enables monitoring of wound healing while being waterproof.



Not made with natural rubber latex.

Features & Benefits:

- Transparent film for easy wound monitoring
- Semipermeable to allow moisture vapor and oxygen exchange
- Includes an absorbent pad for additional fluid management
- Sterile and adhesive for secure placement
- Ideal for cuts, abrasions, and post-operative wounds

Visit Us Online at:
www.ProCureProducts.com

Available Sizes:

Products	Size	Reorder #	Units
Transparent Film with Island Dressing	2.75in x 4in	PCHT124	1/Pouch, 50 Pouches/Box, 40 Boxes/Case
Transparent Film with Island Dressing	2.375in x 2.75in	PCHT122	1/Pouch, 50 Pouches/Box, 40 Boxes/Case
Transparent Film with Island Dressing	3.5in x 4in	PCHT134	1/Pouch, 25 Pouches/Box, 40 Boxes/Case
Transparent Film with Island Dressing	6in x 6in	PCHT166	1/Pouch, 25 Pouches/Box, 40 Boxes/Case