


Transparent Film Roll *Semipermeable Adhesive Film with Grid*



ProCure Transparent Film Roll, is a semipermeable adhesive film supplied in a convenient roll with a printed grid for precise measurement and cutting. Designed to provide a breathable barrier, it allows moisture vapor exchange while protecting against water, bacteria, and contaminants. Flexible and conformable, it adheres smoothly to contoured areas, making it ideal for securing dressings, IV sites, and creating custom-sized coverings for a variety of clinical needs.

**ENHANCED
PATIENT COMFORT**
Featuring a hypoallergenic adhesive film that is gentle on the skin, and allows moisture/air to pass through. Easy-release S-curve liner for hassle-free application.

Not made with natural rubber latex.



Features & Benefits:

- Semipermeable to promote moisture and oxygen exchange
- Transparent for easy monitoring of the wound
- Grid design for easy cutting and customization
- Waterproof and sterile for optimal protection
- Flexible, breathable material for comfort

Visit Us Online at:
www.ProCureProducts.com

Available Sizes:

Products	Size	Reorder #	Units
Transparent Film Roll	2in x 11in	PCHTR211	1 Roll/Box, 90 Boxes/Case
Transparent Film Roll	4in x 11in	PCHTR411	1 Roll/Box, 90 Boxes/Case