

Silver Collagen Matrix Wound Dressing

ProCure Silver Collagen Matrix Wound Dressing is a sterile, advanced wound care dressing composed of collagen, sodium alginate, carboxymethylcellulose (CMC), EDTA, and silver chloride. This specialized matrix supports the body's natural healing process while helping inhibit microbial growth within the dressing. Upon contact with wound exudate, it transforms into a soft gel sheet that maintains a moist environment to promote granulation and epithelialization. The EDTA helps maintain dressing integrity, while the flexible design allows it to be easily trimmed and layered for managing deep or irregularly shaped wounds.

Supports the body's natural wound healing process while gently covering wounds, and maintaining moisture balance for wound closure.



Features & Benefits:

- Sterile, collagen-based matrix with EDTA and antimicrobial silver chloride
- Helps inhibit microbial growth within the dressing
- Forms a soft gel upon contact with wound exudate
- Maintains a moist environment to support tissue regeneration
- Promotes granulation and epithelialization for wound closure
- Can be trimmed and layered for deep or irregular wounds
- Provides comfort, flexibility, and protection throughout healing

DEPTH/STAGE
PARTIAL/FULL
1 2 3 4

ABSORBENCY
MEDIUM

Available Sizes:

Products	Size	Reorder #	Units
Silver Collagen Matrix Dressing	2in x 2in	PCHCMLS22	10/Box, 10 Boxes/Case
Silver Collagen Matrix Dressing	4in x 4in	PCHCMLS44	10/Box, 10 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com

*Not made with
natural rubber latex.*