

Collagen Matrix Wound Dressing

ProCure Collagen Matrix Wound Dressing is a sterile, advanced wound care dressing composed of collagen, sodium alginate, carboxymethylcellulose (CMC), and EDTA. This unique matrix supports the body's natural healing process by maintaining moisture balance and gently covering the wound. Upon contact with wound exudate, it transforms into a soft gel sheet that promotes granulation and epithelialization for optimal wound closure. The dressing can be easily trimmed and layered to manage deep or irregularly shaped wounds, providing flexibility, comfort, and reliable protection throughout the healing process.

Supports the body's natural wound healing process while gently covering wounds, and maintaining moisture balance for wound closure.



Features & Benefits:

- Sterile, collagen-based matrix with EDTA for advanced wound healing
- Forms a soft gel on contact with wound exudate
- Maintains a moist environment to support tissue regeneration
- Promotes granulation and epithelialization for wound closure
- Can be trimmed and layered for deep or irregular wounds
- Helps maintain dressing integrity and comfort during wear



Available Sizes:

Products	Size	Reorder #	Units
Collagen Matrix Dressing	2in x 2in	PCHCMD22	10/Box, 10 Boxes/Case
Collagen Matrix Dressing	4in x 4in	PCHCMD44	10/Box, 10 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com

*Not made with
natural rubber latex.*