

Waterproof Foam Dressing with Adhesive Border

ProCure Waterproof Foam Dressing with Adhesive Border is a moderately absorbent, waterproof sterile foam dressing. It is structured to retain moderate exudate, while maintaining a moist wound environment.

The dressing has a waterproof adhesive backing and a sterile foam cushioning that easily adheres to the skin and seals the wound bed, protecting against bacterial contamination. The waterproof outer layer of the dressing prevents water or other liquids from penetrating the wound site, making it ideal for use in wet environments or for patients who may need to shower or bathe while wearing the dressing. The adhesive border provides a secure and comfortable fit, helping to keep the dressing in place and preventing it from slipping or rolling.

Our Waterproof Bordered Foam Wound Dressing with Adhesive Border is available in a round, fenestrated version with a radial slit and starburst opening for convenient application to tube sites.


ENHANCED PATIENT COMFORT

Multi-layered, waterproof composite dressings provide high absorbency and gently conform to hard-to-dress areas. Designed to protect, cushion, and support optimal healing.

**Special
Size!**

Available Sizes:

Products	Size	Reorder #	Units
Waterproof Foam Dressing	4in x 4in	PCHBF44	10/Box, 40 Boxes/Case
Waterproof Foam Dressing	6in x 6in	PCHBF66	10/Box, 40 Boxes/Case
Waterproof Foam Dressing-Fenestrated	4in Round	PCHFRF4	10/Box, 40 Boxes/Case



Features & Benefits:

- Waterproof barrier protects against fluids and bacteria
- Moderate absorbency helps maintain moisture balance
- Flexible, comfortable fit with secure adhesive border
- Breathable, low-adherent pad for gentle removal
- Available in a round, fenestrated version with radial slit and starburst opening for tube access
- Suitable for showering or bathing

DEPTH/STAGE
PARTIAL/FULL
1 2 3 4

ABSORBENCY
MEDIUM

Visit Us Online at:
www.ProCureProducts.com

*Not made with
natural rubber latex.*