

### Safety LANCETS *BUTTON ACTIVATED*

ProCure Safety Lancets provide a safe, dependable, and easy-to-use option for capillary blood sampling, making them an essential tool in diabetes care. These sterile, single-use lancets feature a retractable needle design that helps prevent accidental needlestick injuries and cross-contamination, ensuring both patient and caregiver safety. Designed with comfort and reliability in mind, they deliver consistent results with every use. The convenient **button-activated style** makes sampling quick and hassle-free, whether in clinical settings or for at-home monitoring.



#### Features & Benefits:

- Available in button-activated style
- Retractable needle design helps prevent needlestick injuries
- Single-use, sterile lancets reduce risk of cross-contamination
- Consistent depth control for reliable blood sampling
- Minimal discomfort for improved patient experience
- Quick, one-step activation for ease of use
- No pre-loading or assembly required
- Compact and portable for home or clinical use
- Not made with natural rubber latex

Visit Us Online at:  
[www.ProCureProducts.com](http://www.ProCureProducts.com)

Product	Size	Reorder #	Units
Safety Lancets-Button	21G x 1.8mm	PCSL21B100	100/box 24 box/case
Safety Lancets-Button	25G x 2mm	PCSL25B100 PCSL25B200	100/box 24 box/case 200/box 24 box/case
Safety Lancets-Button	28G x 1mm	PCSL28B100 PCSL28B200	100/box 24 box/case 200/box 24 box/case
Safety Lancets-Button	30G x 0.7mm	PCSL30B100	100/box 24 box/case