

Sterile Fabric Bandages

ProCure Fabric Bandages provide flexible, comfortable protection for everyday cuts, scrapes, and minor injuries. The sterile, soft fabric conforms easily to the skin for a secure fit, while the absorbent, non-stick pad cushions and protects wounds without sticking. Gentle and latex-free, these bandages enhance comfort and are ideal for daily use.



ENHANCED PATIENT COMFORT

Flexible, comfortable fabric material that won't stick to wounds.



Features & Benefits:

- Flexible fabric conforms for a secure, comfortable fit
- Absorbent, non-stick pad protects without sticking to wounds
- Skin-friendly adhesive leaves no sticky residue
- Sterile and latex-free for safe use on sensitive skin
- Ideal for minor cuts, scrapes, and everyday first aid

Available Sizes:

Products	Size	Reorder #	Units
Fabric Bandages	3/4in x 3in	PCBAF75	1/Pouch, 100/Box, 12 Boxes/Case
Fabric Bandages	1in x 3in	PCBAF1	1/Pouch, 100/Box, 12 Boxes/Case
Fabric Bandages	2in x 4in	PCBAF2	1/Pouch, 50/Box, 24 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com

*Not made with
natural rubber latex.*