

Compression Systems 3 Bandage Kits

ProCure offers two- and four-layer compression bandage systems to help manage edema from venous leg ulcers and related conditions. The **2 Layer Compression System** provides effective graduated compression with a cohesive comfort layer, while the **2 Layer Compression LITE** delivers lighter, slimmer compression with foam padding for easy use in standard or orthotic footwear. The **4 Layer Compression System** combines four wraps with adhesive strips for full compression, comfort, and protection. All systems feature absorbent layers, can be worn up to 7 days (with adequate circulation), and come as convenient all-in-one kits to simplify dressing changes.



2 LAYER COMPRESSION Bandage Kit

Layer 1 – Cohesive Comfort Compression Bandage
4in. x 7.1yds. (10cm x 6.5m) stretched

Layer 2 – Cohesive Compression Bandage
4in. x 8.9yds. (10cm x 8.2m) stretched



2 LAYER COMPRESSION LITE Bandage Kit

Layer 1 – Foam Padding Bandage
4in. x 3.8yds. (10cm x 3.5m)

Layer 2 – Compression Bandage
4in. x 5.5yds. (10cm x 5m) stretched

Convenient
all-in-one kits
simplifies
dressing changes

4 LAYER COMPRESSION Bandage Kit

Layer 1 – Padding Bandage
4in. x 3.8yds. (10cm x 3.5m)

Layer 2 – Conforming Bandage
4in. x 4.9yds. (10cm x 4.5m)

Layer 3 – Compression Bandage
4in. x 9.5yds. (10cm x 8.7m)

Layer 4 – Cohesive Bandage
4in. x 6.5yds. (10cm x 6m)

Adhesive Tape Strips



Features & Benefits:

- Graduated compression supports venous blood flow and reduces swelling
- 2 and 4 layer options suit different compression needs
- Absorbent layers manage wound exudate
- Comfortable, flexible, and easy-to-wear materials
- LITE version fits easily into shoes or orthotics
- 4 layer system includes adhesive strips for secure placement
- All-in-one kits simplify dressing changes

Available Sizes:

Products	Reorder #	Units
2-Layer Compression Bandage Kit	PCHCB002	45 Boxes/Case
2-Layer Compression Bandage Kit LITE	PCHCB002L	45 Boxes/Case
4-Layer Compression Bandage Kit	PCHCB004	64 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com

Not made with
natural rubber latex.