

Non-Sterile Elastic Bandages with Self-Closure

ProCure Elastic Bandages with Self-Closure provide comfortable, dependable support for a wide range of uses. Designed for easy application and secure hold—without clips—these bandages deliver consistent compression to help manage swelling, strains, sprains, and soreness. The soft polyester/cotton blend ensures excellent elasticity, breathability, and comfort. Each bandage is individually poly-wrapped for hygiene and reusability, making it ideal for supporting dressings or securing ice and hot packs.



**Strong
Hold
Reusable**

Body Part	Size of Wrap
Wrist	2 inch
Foot	2 inch
Elbow	3 inch
Ankle	3 inch
Knee	3 inch
Lower Leg	4 inch
Shoulder	4 inch
Chest	6 inch

Features & Benefits:

- Provides consistent, dependable compression for support and swelling control
- Soft polyester/cotton weave ensures elasticity, recoverability, and breathability
- Self-closure keeps the bandage wrapped and secure without clips
- Strong hold stays in place during daily wear
- Reusable and comfortable for extended wear

Available Sizes:

Products	Size	Reorder #	Units
Elastic Bandages with Self Closure	2in x 5yds	PCAB200	10/Box, 5 Boxes/Case
Elastic Bandages with Self Closure	3in x 5yds	PCAB300	10/Box, 5 Boxes/Case
Elastic Bandages with Self Closure	4in x 5yds	PCAB400	10/Box, 5 Boxes/Case
Elastic Bandages with Self Closure	6in x 5yds	PCAB600	10/Box, 5 Boxes/Case

Visit Us Online at:
www.ProCureProducts.com

Not made with
natural rubber latex.