


Silver Alginate Dressing

Contact Layer for Wounds



ProCure Silver Alginate Dressing is a sterile, soft contact layer composed of calcium alginate, sodium alginate, and antimicrobial silver particles. Made from natural seaweed fibers, it is highly absorbent and forms a gel-like surface over the wound, maintaining a moist environment that supports optimal healing while preventing the dressing from sticking to the wound bed. The infused silver provides antimicrobial protection, helping restrain bacterial growth within the dressing for up to 24 hours, and has been shown effective against common pathogens such as *Staphylococcus aureus*, *Escherichia coli*, and *Pseudomonas aeruginosa*.

Flexible and conformable, this non-adherent contact layer enhances patient comfort, minimizes pain during dressing changes, and can be used with most absorbent layers or compression systems.



ENHANCED PATIENT COMFORT

This soft, non-adherent contact layer is infused with silver and gently conforms to delicate wound beds, minimizing pain during dressing changes and supporting a moist healing environment.

Not made with
natural rubber latex.



Features & Benefits:

- Sterile, silver-infused contact layer
- Flexible, non-adherent, and conformable
- Absorbs exudate and maintains a moist environment
- Helps control bacterial growth for up to 24 hours
- Effective against common pathogens
- Compatible with absorbent layers or compression systems
- Enhances patient comfort during dressing changes

Visit Us Online at:
www.ProCureProducts.com

Product	Size	Reorder #	Units
Silver Alginate Dressing	4in x 4in	PCHSCD44	10/Box, 60 Boxes/Case