

Professional Monitoring BLOOD GLUCOSE METER & TEST STRIPS

The TRUE METRIX[®] PRO Blood Glucose Meter and Test Strips system offers a comprehensive solution for professional blood glucose monitoring, ideal for use in clinical settings. Combining the advanced TRUE METRIX[®] PRO Meter with high-precision test strips, this system utilizes Triple Sense Technology[®] to detect, analyze, and correct variables in each sample—delivering fast, accurate, and reliable results. The meter requires no coding, supports alternate site testing, and stores up to 500 readings, while the test strips require only a tiny 0.5 µL sample and provide results in as little as 4 seconds.

Designed for multiple-patient use with enhanced safety features like strip release and audible fill detection, the TRUE METRIX[®] PRO system ensures ease of use, infection control, and confidence in diabetes management.



TRIPLE SENSE TECHNOLOGY[®] automatically detects, analyzes and corrects variables in each blood sample to ensure proven accuracy and confidence in results.*

Event Tagging - The TRUE METRIX[®] PRO System features advanced event tags that combine lifestyle with management.



Features & Benefits:

- No coding required for easy operation
- Fast results in as little as 4 seconds
- Requires only a 0.5 µL blood sample
- Triple Sense Technology[®] ensures high accuracy
- 500-test memory with date and time stamp
- Alternate site testing option for patient comfort
- Single-use, auto-disabling test strips for added safety
- Strip release button and audible fill detection for hygienic, user-friendly handling

Product	Reorder #	Units
Blood Glucose Meter	PCBGM000	1ea/box, 6 boxes/case
Testing Strips	PCBGM050	50ct/box, 24boxes/case
Hi/L0 Control Solution	PCBGMCS13	1ea(hi&lo)/box, 48 boxes/case

Visit Us Online at:
www.ProCureProducts.com