

The ProCure TRUE METRIX® PRO Competitive Advantage

Product Features	TRUE METRIX® PRO Professional Monitoring Blood Glucose System	Assure® Platinum	Assure® Prism
Sample Size	0.5µL	0.5µL	0.5µL
Test Speed	As fast as 4 seconds	7 seconds	5 seconds
Meal Tags	Pre- and Post-meal	None	Post-meal
Event Tags	4	None	None
Strip release	Yes	Yes	Yes
Alternate Site Testing	Yes	No	Yes
Averaging	7/14/30	7/14/30/90	7/14/90
Ketone Test Alert	Yes	No	No
Download Capabilities	Yes TRUEmanager® PRO Diabetes Management Software	No	No
Testing Reminders	Yes (4)	No	Yes (1)
Memory	500 with date and time	500 with date and time	500 with date and time
BG Range (mg/dL)	20-600	20-600	20-600
Hematocrit Range (%)	20-70%	30-55%	20-60%
Blood Sampling	Fresh capillary whole blood from the finger or forearm, or venous blood drawn in sodium heparin tubes	Capillary	Fresh capillary whole blood samples drawn from the fingertips and alternative sites such as the forearm, palm, thigh and calf.


Quality Control Testing

Control Solution	Level 1, Level 3	Level 1, Level 2	Level 1, Level 2
------------------	------------------	------------------	------------------

Please refer to TRUE METRIX PRO Owner's Booklet for safety information and TRUE METRIX PRO Test Strip Instructions for Use for limitations of use. Competitive information on file.

PROCURE Product Description	Reimbursement Code	ProCure Item#
ProCure TRUE METRIX® PRO Professional Monitoring Blood Glucose Meter - 1/ea, 12ea/cs	56151-1476-01	PCBGM000
ProCure TRUE METRIX® PRO Test Strips 50ct - 1/bx, 24bx/cs	56151-1466-04	PCBGM050
ProCure TRUE METRIX® Control Solution Level 1 & Level 3, 3ml/ea - 1/bx, 48 bx/cs	56151-1450-13	PCBGMCS13

For more information, contact your TwinMed sales rep or call 1-877-TwinMed (894-6633). For clinical questions, please call 1-804-ProCure (776-2873).



TRUE METRIX® PRO
Professional Monitoring
Blood Glucose System

www.trividiahealth.com

PROCURE
PRODUCTS

TRUE METRIX is intended for self-monitoring blood glucose only and not for multiple patient use. Only TRUE METRIX PRO is intended for multiple patient use.
NICO-4767 01/26 © 2026 Trividia Health, Inc. TRUE METRIX and TRUEmanager are trademarks of Trividia Health, Inc. All other product names and brand names are trademarks of their respective owners.